
# CytoFLEX™ Flow Cytometer 中文操作手冊

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


1. **Collection.** Establishes control over data recording options and displays the acquisition status.
2. **Test tubes.** Allows you to configure and duplicate sample tubes, set display attributes, manage experimental data and compensation.
3. **Plot area.** Includes plot and gating controls, as well as an area for drawing plots and generating graphs.

### 機器操控區(Collection)


1. **Acquisition control.** Controls sample loading/unloading and data acquisition and recording.
2. **Acquisition status.** Displays such information as the acquisition rate (Events/Sec), cell count, duration and abort (%).
3. **Acquisition conditions.** Sets the necessary conditions for recording data.
  - **Events to Record.** Used to set the number of events to record in the specified population.
  - **Time to Record.** Used to set the collection time duration in seconds.
4. **Sample flow rate.** Sets the acquisition rate for data collection.
  - **Slow** : 10  $\mu$ L/min    - **Medium** : 30  $\mu$ L/min    - **High** : 60  $\mu$ L/min
  - **Custom** : 10 - 240  $\mu$ L/min

 Initialize	Put the instrument in initialized state.
 Standby	Put the instrument in standby state.
 Run	Start acquisition or continue an acquisition if previously stop.
 Stop	Stop the acquisition of the current sample and output the results.
 Record	Used to set the collection conditions for sample recording.
 Restart	Reset the current acquired events to zero and clear the current data in memory. Acquisition restarts at zero events.
 Backflush	Flush the sample line and flow cell with sheath fluid to remove bubbles.
 Boost	To transfer the sample to the flow cell.
 Next Tube	Switch to the next sample tube.
 Acq. Setting...	Display the acquisition setting dialog box to adjust the Cytometer settings.

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



	New Tube	Create a new tube.
	Duplicate without data	Create a blank row at the end of the current experiment.
	Delete Tube	Delete one currently selected tube.
	Delete Multiple Tubes	Delete the multiple selected tubes.
	Property	Sample tube basic information.
	Open the folder	Open the experiment folder.
	Compensation Matrix	Open the compensation matrix.

白	Indicates that the tube data was not collected.
藍	Indicates that the tube data was acquired by selecting <b>Run</b> but not saved, and the tube data can be overwritten.
綠	Indicates that the tube data was saved by selecting <b>Record</b> and that this data cannot be overwritten.

## 繪圖區(Plot area)

	Histogram	Create a Histogram Plot and specify the plot properties.
	Dot Plot	Create a Dot Plot and specify the plot properties.
	Density Plot	Create a Density Plot and specify the plot properties.
	Pseudo Color Plot	Create a Pseudo Color Plot and specify the plot properties.
	Contour Plot	Create a Contour Plot and specify the plot properties.
	New Histograms	Create multiple Histogram Plots and specify the plot properties.
	New 2-D Plots	Create multiple Dot Plots and specify the plot properties.
	Statistics	Create Statistical charts.
	Population Hierarchy	Create Hierarchical charts.
	Line Segment	Insert a Linear gating of plots.
	Vertical	Insert a Vertical gating of plots.
	Lasso	Insert a Lasso gating into a dual parameter plots.
	Polygon	Insert a Polygon gating into a dual parameter plots.
	Rectangle	Insert a Rectangle gating into a dual parameter plots.
	Four Quadrant	Insert a Four Quadrant gating into a dual parameter plots.
	Hinged	Insert a Hinged gating into a dual parameter plots.
	Auto Line Segment	Crear an Auto Line Segment around the selected population on a plot.
	Auto Polygon	Crear an Auto Polygon around the selected population on a plot.
	Zoom In	For Zooming in.
	Zoom Out	For Zooming out.
	Pan	For scaling axis ranges in the plots.
	Single Side Pan	For scaling single axis range in the plots.
	Adjust Gain	For increasing and lowering gain adjustments on the plots.
	Adjust Compensation	For adjusting compensation of either of the parameters on a 2-D histogram.
	Threshold	For setting the minimum particle size limit or fluorescence intensity that acquisition will allow.
	Undo	For undoing an action in the drawing area.
	Redo	For redoing an action in the drawing area.

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- The gain differences must be  $\leq 20\%$  from the target gain.
- The median fluorescence intensity (MFI) differences must be  $\leq 5\%$  from the target MFI.
- The rCV must be  $\leq 5\%$ .

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



4. **Collection.** Establishes control over data recording options and displays the acquisition status.
5. **Test tubes.** Allows you to configure and duplicate sample tubes, set display attributes, manage experimental data and compensation.
6. **Plot area.** Includes plot and gating controls, as well as an area for drawing plots and generating graphs.

### 機器操控區(Collection)


1. **Acquisition control.** Controls sample loading/unloading and data acquisition and recording.

- 2. Acquisition status.** Displays such information as the acquisition rate (Events/Sec), cell count, duration and abort (%).
- 3. Acquisition conditions.** Sets the necessary conditions for recording data.
- **Events to Record.** Used to set the number of events to record in the specified population.
  - **Time to Record.** Used to set the collection time duration in seconds.
- 4. Sample flow rate.** Sets the acquisition rate for data collection.
- **Slow** : 10  $\mu\text{L}/\text{min}$     - **Medium** : 30  $\mu\text{L}/\text{min}$     - **High** : 60  $\mu\text{L}/\text{min}$
  - **Custom** : 10 - 240  $\mu\text{L}/\text{min}$

 Initialize	Put the instrument in initialized state.
 Standby	Put the instrument in standby state.
 Run	Start acquisition or continue an acquisition if previously stop.
 Stop	Stop the acquisition of the current sample and output the results.
 Record	Used to set the collection conditions for sample recording.
 Restart	Reset the current acquired events to zero and clear the current data in memory. Acquisition restarts at zero events.
 Backflush	Flush the sample line and flow cell with sheath fluid to remove bubbles.
 Boost	To transfer the sample to the flow cell.
 Next Tube	Switch to the next sample tube.
 Acq. Setting...	Display the acquisition setting dialog box to adjust the Cytometer settings.
 Eject	Open the plate holder.
 Auto Record	Sample acquisition occurs in the order indicated by the numbers

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



1. **Test tube management controls.** Manages sample tubes. Used to add, copy, or delete attributes, open the tube property, and open the compensation matrix.
2. **Test tube status indication.** Displays a colored symbol in front of each tube indicating the status of the tube processing.
3. **Test tube list.** Displays the sample tubes used in the experiment. Right-click a tube in the list to perform additional operations.

	Plate	Open the plate window.
	Add Plate	Add one new plate.
	Add Plate from Template	Add a plate template with preset settings.
	Duplicate Current Plate without Data	Creat a copy of the selected plate without data.
	Delete Plate	Delete one plate.
	Save Plate Template	Save the plate condition as a template.
	Set As Sample Wells	Set selected wells as sample wells.
	Set As Cleaning Agent Wells	Set selected wells as cleaning wells.
	Set As Deionized Water Wells	Set selected wells as Deionized Water wells.
	Set As Empty Wells	Reset selected wells as empty.
	Set Acquisition Condition	Select the desired acquisition settings.
	Set Auto Acquisition	Set the selected wells for auto record.
	Cancel Auto Acquisition	Remove the auto record setting from the selected wells.
	Heat Map	Open Heat Map window.
	New Heat Map	Add one new plate with Heat Map.


	Duplicate	Duplicate current Heat Map.
	Modify Current Heat Map Setting	Modify existing Heat Map settings.
	Delete Current Heat Map	Delet a single Heat Map from the list of Heat Map window.
	Delete Multiple Heat Maps	Delet multiple Heat Maps from the list of Heat Map window.
	Refresh	Refresh a single Heat Map from the list of Heat Map window.
	Refresh All	Refresh all Heat Maps from the list of Heat Map window.
	Export to Graphic File	Export a Heat Map as a Graphics File(.bmp or .emf)
	Export to Clipboard	Export a Heat Map to a Clipboard file(.bmp)

 藍勾	A blue check mark means the data is acquired, not recorded.
 綠勾	A green check mark means the data is recorded.
 紅叉	A red cross mark means the data is recorded. But the acquisition is terminated abnormally. For example, the well is skipped or the acquisition is manually stopped.

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Acq. Setting：調整 Gain、Threshold 及 Width。

Set As Sample Wells

Naming Rules  
Prefix: Tube  
Example: 01-Tube-A1

Mix/Backflush  
 Mix 3 sec (0.1 ~ 100)  
 Backflush 3 sec (0.1 ~ 100)

Group  
Name:   Apply to Sample ID  
Color:

Acq.Setting Channel Compensation Matrix Stopping Rules


Gain  
FSC 500 (1~3000)  
SSC 500 (1~3000)  
FITC 500 (1~3000)  
PE 500 (1~3000)  
ECD 500 (1~3000)  
PC5.5 500 (1~3000)  
PC7 500 (1~3000)  
APC 500 (1~3000)  
APC-A700 500 (1~3000)  
Default Recommended

Threshold  
Primary Threshold (Trigger Level)  
Channel: FSC  
 Manual 10000 (>0)  Height  Area  
 Automatic  
Logic Operator:   
Secondary Threshold (Trigger Level)  
Channel:   
 Manual 1 (>0)  Height  Area  
 Automatic  
Default

Width  
Channel: FSC  
Default

Import from File... Import from Catalog...  
Import Acquisition Conditions from FCS File...  
OK Cancel

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Set As Sample Wells

Naming Rules  
Prefix: Tube  
Example: 01-Tube-A1

Mix/Backflush  
 Mix 3 sec (0.1 ~ 100)  
 Backflush 3 sec (0.1 ~ 100)

Group  
Name:   Apply to Sample ID  
Color:

Acq.Setting Channel Compensation Matrix Stopping Rules

Use	Channel	Label
<input checked="" type="checkbox"/>	FITC	
<input checked="" type="checkbox"/>	PE	
<input type="checkbox"/>	ECD	
<input checked="" type="checkbox"/>	PC5.5	
<input checked="" type="checkbox"/>	PC7	
<input checked="" type="checkbox"/>	APC	
<input checked="" type="checkbox"/>	APC-A700	
<input checked="" type="checkbox"/>	APC-A750	
<input checked="" type="checkbox"/>	PB450	
<input checked="" type="checkbox"/>	KO525	
<input checked="" type="checkbox"/>	Violet610	
<input type="checkbox"/>	Violet660	
<input type="checkbox"/>	Violet780	

Import Acquisition Conditions from FCS File...  
OK Cancel

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Mix/Backflush  
 Mix 3 sec (0.1 ~ 100)  
 Backflush 3 sec (0.1 ~ 100)

Group  
Name:   
 Apply to Sample ID  
Color:

Acq.Setting Channel Compensation Matrix Stopping Rules

Use  Show Autofluorescence  Area and Height in Sync

Autofl.	Channel	-FITC%	-PE%	-PC5.5%	-PC7%	-APC%	-APC-A700%	-APC-A750%	-PB450%	-KO525%	-Violet610%
0.00	FITC		0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00
0.00	PE	0.00		0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00
0.00	PC5.5	0.00	0.00		0.00	0.00	0.00	0.00	0.00	0.00	0.00
0.00	PC7	0.00	0.00	0.00		0.00	0.00	0.00	0.00	0.00	0.00
0.00	APC	0.00	0.00	0.00	0.00		0.00	0.00	0.00	0.00	0.00
0.00	APC-A700	0.00	0.00	0.00	0.00	0.00		0.00	0.00	0.00	0.00
0.00	APC-A750	0.00	0.00	0.00	0.00	0.00	0.00		0.00	0.00	0.00
0.00	PB450	0.00	0.00	0.00	0.00	0.00	0.00	0.00		0.00	0.00
0.00	KO525	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00		0.00
0.00	Violet610	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00	

Import from Library... Import...

Import Acquisition Conditions from FCS File...

OK Cancel

### 4. Stopping Rules : 設定儲存條件。

Set As Sample Wells

Naming Rules  
Prefix:   
Example: 01-Tube-A1

Mix/Backflush  
 Mix 3 sec (0.1 ~ 100)  
 Backflush 3 sec (0.1 ~ 100)

Group  
Name:   
 Apply to Sample ID  
Color:

Acq.Setting Channel Compensation Matrix Stopping Rules

Events to Record:  Events  
in

Time to Record:  sec



Import Acquisition Conditions from FCS File...

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
OK Cancel

## 繪圖區(Plot area)

	Histogram	Create a Histogram Plot and specify the plot properties.
	Dot Plot	Create a Dot Plot and specify the plot properties.
	Density Plot	Create a Density Plot and specify the plot properties.
	Pseudo Color Plot	Create a Pseudo Color Plot and specify the plot properties.
	Contour Plot	Create a Contour Plot and specify the plot properties.
	New Histograms	Create multiple Histogram Plots and specify the plot properties.
	New 2-D Plots	Create multiple Dot Plots and specify the plot properties.
	Statistics	Create Statistical charts.
	Population Hierarchy	Create Hierarchical charts.
	Line Segment	Insert a Linear gating of plots.
	Vertical	Insert a Vertical gating of plots.
	Lasso	Insert a Lasso gating into a dual parameter plots.
	Polygon	Insert a Polygon gating into a dual parameter plots.
	Rectangle	Insert a Rectangle gating into a dual parameter plots.
	Four Quadrant	Insert a Four Quadrant gating into a dual parameter plots.
	Hinged	Insert a Hinged gating into a dual parameter plots.
	Auto Line Segment	Creat an Auto Line Segment around the selected population on a plot.
	Auto Polygon	Creat an Auto Polygon around the selected population on a plot.
	Zoom In	For Zooming in.
	Zoom Out	For Zooming out.
	Pan	For scaling axis ranges in the plots.
	Single Side Pan	For scaling single axis range in the plots.
	Adjust Gain	For increasing and lowering gain adjustments on the plots.
	Adjust Compensation	For adjusting compensation of either of the parameters on a 2-D histogram.
	Threshold	For setting the minimum particle size limit or fluorescence intensity that acquisition will allow.
	Undo	For undoing an action in the drawing area.
	Redo	For redoing an action in the drawing area.

	Align Left	Align all the selected items to the left of the selection area.
	Align Right	Align all the selected items to the right of the selection area.
	Align Top	Align all the selected items to the top of the selection area.
	Align Bottom	Align all the selected items to the bottom of the selection area.
	Vertical Distribute	Align all the selected items to the vertical distribution.
	Horizontal Distribute	Align all the selected items to the horizontal distribution.
	Make Same Width	Resize the selected items to all be the same width as the reference item.
	Make Same Height	Resize the selected items to all be the same height as the reference item.
	Make Same Size	Resize the selected items to all be the same size as the reference item.
	Rearrange	For restoring the plots to the default positions.
	Print	For printing and previewing the plot area.
	Print Preview	Used to access the Preview screen.
	Page Setup	Used to adjust the page settings.
	Batch Print	Used to print data for multiple tubes.
	Batch Export to PDF File	Used to print a PDF of the data for multiple tubes.

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The diagram illustrates the transformation of a CSV file. On the left, the 'Old' CSV content is shown as a list of rows with columns: WellLabel, TubeName, SampleID, Group, Name1, Name2, Name3. Below it is an Excel spreadsheet with columns A-H and rows 1-12. On the right, the 'New' CSV content is shown with an additional column Name4 and updated data for rows 1-5, plus a new row 6. Below it is an Excel spreadsheet with columns A-I and rows 1-13, where row 6 is highlighted in green. A purple arrow points from the 'Old' state to the 'New' state.

b. 點擊Add Plate from Layout...，選擇已修改過的001.csv檔並點選Open。

The image shows two screenshots from a software interface. The top screenshot shows a 'Plate' configuration window with a dropdown menu open, highlighting the 'Add Plate from Layout...' option. The bottom screenshot shows a Windows file explorer window with the file '001.csv' selected in the 'Test' folder. The file name '001.csv' and file type 'CSV File (\*.csv)' are visible at the bottom of the window.

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#### Materials :

1. Propidium iodide (sigma) : 10 mg/ml in water, store at 4 degrees.

Note : solution made at RT will fall out of solution in fridge. No problem. Just mix it up and squirt the cloud into the staining mix. It will go into solution there.

2. RNaseA (Sigma), 10 mg/ml in water, store at 4 degrees (long term 20 degrees)
3. PBS : at 4 degrees
4. Ethanol (100%, store at -20 degrees)

#### Fixation :

- (1) Spin cells out of media – 1500 rpm x 5 minutes

- (2) Wash once with PBS
- (3) Collect  $2 \times 10^6$  cells.
- (4) Pellet cells by spinning at 1500 rpm, 4°C for 5 minutes.
- (5) Resuspend cell pellet in 300ul of cold PBS. Vortex.
- (6) Fix cells by adding 700ul of -20°C absolute ethanol. (drop by drop initially)
- (7) Store cells at -20°C in this fixation buffer until ready for analysis. (No more than 2 weeks)

Staining :

- (1) Centrifuge (as above) fixed cells and resuspend pellet in 1 ml of PBS.
- (2) Add 100ul of 200 ug/ml DNase-free, RNase-A and incubate at 37°C for 30 minutes.
- (3) Add 100ul of 1 mg/ml propidium iodide (light sensitive) and incubate at room temperature for 5-10 minutes.
- (4) Place samples in 12x75mm tubes and read on Flow cytometry

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